

Product Comparison

AED Selection Guide



LIFEPAK® 1000 Defibrillator



LIFEPAK CR® Plus Defibrillator

	Manual Mode	AED Mode	AED Mode
PRIMARY USE	PROFESSIONAL BLS/ ALS RESPONDERS	BLS FIRST RESPONDERS	LAY OR TRAINED RESPONDERS
Device users	Medical Clinics, EMS	EMS, Fire Service, Police, Surf Lifesaving, Industrial or frequent users (> yearly)	Commercial or infrequent users (< yearly)
Energy Sequence – AED	User configurable to meet ARC/NZRC guidelines 200J – max energy 360J	User configurable to meet ARC/NZRC guidelines 200J – max energy 360J	Configured to meet ARC/NZRC guidelines 200J – max energy 360J
THERAPY FEATURES			
Defibrillation	Escalating Energy 200 – 360J	Escalating Energy 200 – 360J	Escalating Energy 200 – 360J
Manual Shock	Yes	N/A	N/A
Semi-automatic AED Shock	Yes	Yes	Yes
MONITORING FEATURES			
ECG Monitoring	Yes	N/A	N/A
DEVICE FEATURES			
Voice Prompts	Yes	Yes	Yes
Readiness Indicator	OK, BATTERY, SERVICE	OK, BATTERY, SERVICE	OK, BATTERY, ATTENTION, SERVICE
Electrode Options	EDGE System™ pre-connected electrodes Paediatric (Infant/Child) Reduced Energy Defibrillation Electrodes	EDGE System pre-connected electrodes Paediatric (Infant/Child) Reduced Energy Defibrillation Electrodes	QUIK-PAK™ electrodes with REDI-PAK™ preconnect system Paediatric (Infant/Child) Reduced Energy Defibrillation Electrodes
ECG Cables	3 wire	N/A	N/A
Lead Views	II	N/A	N/A
Dust/Water Rating	IP55 (Protection from jetting water)	IP55 (Protection from jetting water)	IPX4 (protection from splashing water)
Post Event Review	Infrared wireless to CODE-STAT™ software	Infrared wireless to CODE-STAT software	Infrared wireless to CODE-STAT software
Power Type	Battery – 17 hours operating time or 440 shocks at 200J	Battery – 17 hours operating time or 440 shocks at 200J	SafeGuard™ Power System Rechargeable Battery – 30 discharges or 210 minutes of on time
Included Accessories	Carry Case, Spare electrodes, ECG cable case	Carry Case, Spare electrodes	Carry Case, Spare electrodes, Ambu® Res-Cue Mask® Kit
Warranty	5 Years	5 Years	8 Years

For further information please contact your local Physio-Control representative or visit our website at www.physio-control.com



Physio-Control Headquarters
Redmond, WA 98052 USA
Tel 425 867 4000
Fax 425 867 4121
www.physio-control.com

**Physio-Control
Australia Pty Ltd**
Suite 4.01
15 Orion Road
Lane Cove
NSW 2066
Australia
Toll Free Tel 1800 987 982
Toll Free Fax 1800 890 892



Physio-Control, Inc., 11811 Willows Road NE, Redmond, WA 98052 USA
